
BabyFix Cocoon

Instructions for Use
Gebrauchsanweisung
Mode d'emploi
Istruzioni per l'uso
使用説明書
Instrucciones de uso
Instruções de utilização
Brugsanvisning
Käyttöohje
Bruksanvisning
Gebruiksaanwijzing
Инструкция по эксплуатации
Instrukcja obsługi
사용 설명서
使用説明
Kullanım Kitabı
Návod k použití


Positioning aid for the immobilization of newborns (2 to 5 kg) during MR imaging

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Short-term patient positioning and stabilization device for Neonates (2-5kg) during invasive and non-invasive procedures to be used on intact patient's skin. Non-Sterile and not to be used in the sterile field.

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
1 Spread open the BabyFix Cocoon. Make sure that the pearls are evenly distributed. Place the dressed newborn or cotton towel wrapped newborn in the positioning aid. Make sure the shoulders are flush with the mattress and the head is centered.

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2 Close the BabyFix Cocoon keeping it tight fitting. When closing the BabyFix Cocoon make sure the clothing or the cotton towel is as wrinkle-free as possible to avoid pressure points. For sedated newborns the tubes and cables are secured by the tube stop. Use the viewing window to check respiratory movements and attached monitoring systems.

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
3 Place ear protection as usual. Secure the head section of the positioning aid with a fleece strip over the forehead of the newborn.

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4 Connect the vacuum hand pump to the tube and extract as much air as possible from the positioning aid. Close the hose clamp and remove the pump. Tighten the fleece strip on the head if necessary.

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5 Position the newborn in the BabyFix Cocoon in the MR Scanner.

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6 After the exam has been carried out (observe the instructions of the appliance manufacturer) unfasten the fleece strip and let air flow into the positioning aid again by opening the hose clamp.

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
Clean the positioning aid with conventional disinfectants. Apply the disinfectant according to the manufacturer's instructions and wipe clean.

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
Note: Do not immerse product in cleaning fluid.

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Contraindications: Do not use the product for purposes other than those described.

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MR MR safe

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NON STERILE Do not sterilize. A hygiene barrier (mechanical barrier) must always be used. Do not use in case of open wounds. Avoid contact with sharp objects. Remove the pump from the positioning aid during the exam. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare professional (Rx only).

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