

BD Connecta™ Stopcock with extension tube
Sterile
394995 – 394997 – 394926 – 394951 – 394961 –
394982 - 394983

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins, Switzerland
bd.com

TDS number: V201-033 – Rev.02
2020-May

1. General Information

1.1 Intended use


BD Connecta™ Stopcock with extension tube is designed to supply the blood system with fluids or medication coming from one or two different sources via an IV cannula or extension tube. The product can also be used in different combinations in several applications of nutrition administration.

BD Connecta™ Stopcock with Low Volume Extension tube (for pediatric use) is designed to supply the blood system with fluids or medication from one or two different sources via an IV cannula or an extension tube. Alternatively, one of the connections can be used for injections or sampling. The product can also be used in various combinations for several nutrition administration applications.

BD Connecta™ Stopcock with Low Volume Extension is suitable for applications where a stopcock is preferably placed at some distance from the insertion site.

1.2 General description

BD Connecta™ Stopcock with extension tube has a transparent housing with a white or blue tap. The housing has two female Luer-lock fittings with white plugs, and one fitting to which an extension line or a low volume extension tube is assembled. The luer cone at the end of the extension line is protected by a white protection cap. The rotating tap has closed and open position indicators when snap turned (ON-OFF). The rotating tap shows which ports are open. The tap rotation is 360°.



BD Catalog Number	BD Product Description	Tube Length (cm)	Tap color
394995	CONNECTA PLUS3 10CM WHITE	10	White
394997	CONNECTA PLUS3 10CM BLUE	10	Blue
394926	CONNECTA WHT 360DEG TB 25CM	25	White
394951	CONNECTA WHT 360DEG TB 50CM	50	White
394961	CONNECTA WHT 360DEG TB 100 CM	100	White
394982	CONNECTA WHT 360DEG LOW VOLUME TB 15CM	15	White
394983	CONNECTA WHT 360DEG LOW VOLUME TB 30CM	30	White

Please check BD catalog number availability in your country.

Further features: N/A

1.3 Certification

BD Catalog Number	BD Legal Manufacturer and ISO 13485 Certification	CE Certificate Number And Notified Body Brief Name	BD Manufacturing Site (Country of Origin) and ISO 13485 Certification	EC Representative (if applicable)
394995 394997 394926 394951 394961 394982 394983	Address: Becton Dickinson Infusion Therapy AB Florettgatan 29C PO Box 631 SE-251 06 Helsingborg Sweden ISO 13485 Certificate No.: MD 597883	CE certified with BSI (2797) Certificate No.: CE 597884	Address: Becton Dickinson Infusion Therapy systems Inc. S.A. de C.V. Periferico Luis Donaldo Colosio #579 Nogales, Sonora C.P. 84048 Mexico ISO 13485 Certificate No: FM 673986	N/A

1.4 Materials

Component	Material
Housing	Polycarbonate
Tap	Polyethylene
Tube	Polyvinylchloride
Nut	Polycarbonate
Plug	Polypropylene
Protection cap	Polyethylene Master batch white
Silicone liquid	Polydimethylsiloxane
Peg	Polyethylene

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization.
The content of this document may be subject to change without notice.

1.5 Materials of concern

Materials of concern are chemicals or substances that have been identified as having the potential to cause long term effects on humans or the environment.

Material	Comment
Phthalates	<p>Based on our ongoing data collection efforts and/or information received from our suppliers as per 13 March 2019, BD has not identified any</p> <ul style="list-style-type: none"> • Bis(2-ethylhexyl) phthalate (DEHP) (CAS# 117-81-7), • Dibutyl phthalate (DBP) (CAS# 84-74-2), • Benzyl butyl phthalate (BBP) (CAS# 85-68-7), • Bis(2-methoxyethyl) phthalate (DMEP) (CAS# 117-82-8), • Diisobutyl phthalate (DIBP) (CAS# 84-69-5), • Diisopentyl phthalate (DIPP) (CAS# 605-50-5), • Dipentyl phthalate (DPP) (CAS# 131-18-0), • Di-n-hexyl phthalate (DnHP) (CAS# 84-75-3), or • N-pentyl-isopentylphthalate (CAS# 776297-69-9) <p>in the articles and packaging with the product numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w). Furthermore, these products are not designed to contain any phthalates.</p>
Latex	<p>Based on our ongoing data collection efforts and/or information received from our suppliers as per 02 February 2019, the articles with the product numbers above are not formulated with natural rubber latex.</p>
Bisphenol A	<p>Based on our ongoing data collection efforts and/or information received from our suppliers as per 13 March 2019, BD has not identified any</p> <ul style="list-style-type: none"> • 4,4'-isopropylidenediphenol (BPA) (CAS# 80-05-7) <p>in the articles and packaging with the product numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w).</p> <p>There are polycarbonate components in these products. Bisphenol A (BPA), CAS# 80-05-7, is an organic compound that is a chemical building block for polycarbonate. Based on information from our suppliers and BD test results, the BPA level is less than 0.1% wt/wt <1000 ppm). These levels are below a de minimis concentration with no demonstrable clinically significant exposure nor toxicity.</p> <p>No labeling for California Prop 65 is needed.</p> <p>No REACH SVHC declaration is required.</p>
Substances of animal origin BSE/TSE	<p>The raw materials used in the manufacture of these devices do not contain any animal tissue but may contain very small amounts of chemicals that are derived from animal-derived raw materials. These products are manufactured using polymer resins which may contain very small amounts of stearic acids and related substances derived from tallow derivatives. Our resin suppliers have confirmed that these tallow-derived materials have been produced with multiple cycles of conditions at least as rigorous as those specified in Annex C.5 of EN 22442-1 and Section 6 of EMA 410/01 Rev. 3.</p> <p>Therefore, these raw materials meet or exceed the requirements of EN ISO 22442-1 and EMA 410/01 Rev. 3. Based on this information, these products are considered not to present any risk with respect to TSE/BSE or other animal-borne diseases.</p> <p>Furthermore, as recognized by MEDDEV 2.4/1, such derivative chemicals produced in accordance with the aforementioned standards and guidelines are considered irrelevant when determining the classification of a medical device (per MDD 93/42/EEC and EU No 722/2012).</p>
Polyvinyl chloride (PVC)	<p>The medical devices referenced above do contain PVC which is not formulated with any regulated phthalates.</p>

1.6 **REACH information**

Based on our ongoing data collection efforts and/or information received from our suppliers as per 13 March 2019, BD has not identified any chemicals in the articles and packaging of BD Connecta™ Stopcock with extension, in an individual concentration above 0.1% weight by weight (w/w), which have been listed as SVHC and included in the "Candidate list" published by the European Chemical Agency (ECHA) on 27 June 2018 according to Art. 59 (1,10) of the Regulation (EC) N°1907/2006 (REACH).

1.7 **Biocompatibility**

BD Medical products comply with the requirements of the standard for toxicity, pyrogenicity and biocompatibility of medical devices, ISO 10993 series - Biological Evaluation of Medical Devices.

1.8 **Sterilization method**

Sterilization method: **Radiation** (EN ISO 11137-1 "Sterilization for Healthcare products- Radiation –Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices).

1.9 **Shelf life and storage conditions**

The BD Connecta™ Stopcock with extension tube shelf life has been assessed by stability studies in order to verify the functionality, physico-chemical and microbial properties over time.

BD Connecta™ Stopcock with extension tube have a shelf life of 3 years.

Store in a dry and warm place, not exposed to strong light.

1.10 **Standards**

As per extract from the Declaration of Conformity:

Harmonized Standards	
EN 556-1:2001	Sterilization of Medical Devices - Requirements for Medical Devices to Be Designated "Sterile" - Part 1: Requirements for Terminally Sterilized Medical Devices
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN ISO 13485:2016	Medical devices - Quality management systems - Requirements for regulatory purposes.
EN ISO 14971:2012	Medical devices. Application of risk management to medical devices.
EN 20594-1:1993	Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment.
EN ISO 10993-1:2009	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2:2006	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization.

The content of this document may be subject to change without notice.

Harmonized Standards	
EN ISO 11137-1:2006	Sterilization of health care products -- Radiation -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 11137-2:2013	Sterilization of health care products -- Radiation -- Part 2: Establishing the sterilization dose
EN ISO 11737-1:2006	Sterilization of medical devices – Microbial methods- Part 1 : Determination of a population of microorganisms on products
EN ISO 11737-2:2009	Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the validation of a sterilization process
EN 980:2008	Symbols for use in the labelling of medical devices
EN 62366:2008	Medical devices. Application of usability engineering to medical devices
Non-Harmonized Standards	
ISO 14644-1:1999	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
ISO 594-2:1998	Conical fittings with 6% (Luer) taper for syringes, needles and certain other medical equipment –Part 2: Lock Fittings
EN ISO 15223-1:2012	Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General Requirements

Note:

The above standards reflect the status at the time of drafting this document. More information or updates are available on request in the Declaration of Conformity.

1.11 Classification

The products have been classified in accordance with the Medical Devices Directive Annex IX as amended. Rule 2 is applicable and the products belong to Class IIa.

1.12 GMDN code

According to ISO 15225 (Medical devices - Quality management - Medical device nomenclature data structure), BD Connecta™ Stopcock with extension are referenced as follows:

GMDN Code: 32172

GMDN Term: Infusion Stopcock

1.13 Manufacturing practices

The entire manufacturing and testing processes are following the Manufacturing Practices as specified below:

- Incoming raw materials are verified via material inspection and testing and our suppliers are approved via our vendor management system.
- In addition to the automatic on-line inspections, in-process inspections are performed in addition to final product testing to ensure compliance with approved specifications.
- The manufacturing and testing details of each batch of product are recorded on a batch record which is retained in accordance with our document control procedures.
- BD operates a system of internal and external audits to maintain compliance.
- BD confirms that it will continue to adhere to relevant international standards in designing and manufacturing its products.

1.14 Other information

- (Material) Safety Data Sheets are not required for this product.
- Certificate of Food Contact (*Commission Regulation EU 1183/2012 on "plastic materials and articles intended for contact with food" and Directive 2002/72/CE (as amended) "relating to plastic materials and articles intended to come into contact with foodstuffs"*) is not required as BD products are used for general purpose injection and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.
- Good Manufacturing Practices as defined by the FDA Pharmaceutical is not applicable for Medical Devices.

2. Packaging

2.1 Packaging configuration

BD Catalog Number	BD Product Description	Primary Packaging (Qty)	Shelf Box (Qty)	Shipping Case (Qty)	IFU Insert N/A / Yes / No*
394995	CONNECTA PLUS3 10CM WHITE	1	50	250	Yes
394997	CONNECTA PLUS3 10CM BLUE	1	50	250	Yes
394926	CONNECTA WHT 360DEG TB 25CM	1	20	100	Yes
394951	CONNECTA WHT 360DEG TB 50CM	1	20	100	Yes
394961	CONNECTA WHT 360DEG TB 100 CM	1	20	100	Yes
394982	CONNECTA WHT 360DEG LOW VOLUME TB 15CM	1	50	250	Yes
394983	CONNECTA WHT 360DEG LOW VOLUME TB 30CM	1	50	250	Yes

*"No": IFU may be available but not as an insert.


2.2 Packaging material

Component	Material
Unit Pack	Plastic film and medical grade paper
Shelf Box	Cardboard
Shipping Case	Cardboard
IFU	Paper


2.3 Examples of labeling

Labels: According to European Medical Device directive, labels are multilingual.

Primary Packaging Label (Top Web) extracted from document H0570 related to reference 394995:



Shelf Box extracted from document F1922 related to reference 394995:



Shelf label extracted from document D16662 related to reference 394995:

Ref 394995

BD Connecta™

10 cm  **360°**

50

BD Luer-Lok™
Sterile R

Product of Mexico
Сделано в Мексике
Виготовлено в Мексикі

PG CP

11076/2011
16.08.2013


РУ № ФСЗ 2008/03250 от 31.12.2010


(17)YYMM00(10)YJJSSS(30)50


(01)30382903949954

 YYYY-MM  YYYY-MM **LOT** YJJSSS


Shipping Case extracted from document F3742 related to reference 394995:



Case Label extracted from document D16663 related to reference 394995:

Ref 394995

BD Connecta™



BD Luer-Lok™
Sterile[R]

Product of Mexico
Сделано в Мексике
Виготовлено в Мексиці

250



CP	11076/2011
	16.08.2013

РУ № ФСЗ 2008/03250 от 31.12.2010



((17)YYMM00(10)YJJSSS(30)0250



(01)50382903949958

⌚ YYYY-MM 📄 YYYY-MM [LOT] YJJSSS

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization.
The content of this document may be subject to change without notice.

