




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com / en /](http://www.medline.com/en/)



## OPS™ Essential

## Fenestrated drape sheet

### Reference(s)

**ES15230CE**

### Description

A peel pouch, 23x33 cm, contains:  
- 1 bi-laminate drape sheet 45x75 cm, 7 cm diameter fenestration

### Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

### Pack factor

180 units per carton  
60 units per dispenser box

### Packaging

Paper and film peel pouch  
3 dispenser boxes contain all the units in a carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE EO



© 2014 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved..

LIT-ES15230CE-EN01

01.01.INT.12