

**Medline Industries Ltd**

3rd Floor, Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical sheet

High-Absorbency SAP Bed Underpad

Reference(s)

Product No.	Size (cm)	Weight (gm)	Packaging	Min. average absorption (ml) (ISO 11948-1)
MSC282046	58.4 x 58.4	53	100 each (4 bags of 25)	750

Description

- Impervious polypropylene backsheet
- Soft nonwoven topsheet
- All four sides sealed to prevent leakage
- Absorbent layer of cellulose fibers and super absorbent polymer

Main materials

Polypropylene backing
Cellulose
Super absorbent polymer (polyacrylate)

Pack factor

Refer to description item

Legal Manufacturer

Medline Industries LP,
Three Lakes Drive, Northfield Illinois, 60093 USA
EN ISO 13485

Regulation

In compliance with Medical Devices Regulation 2017/745 and UK MD
MD Class: I Non-sterile

- Please refer to complete labeling and instruction for use (if applicable) before using the product.