




Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business
Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

OPS™ Ultimate

Neonatal drape

Reference(s)

TB29494CE

Description

A sterile wrap 60x60 cm contains:
- 1 tri-laminate Neonatal drape, 112x152 cm, 47x51 cm reinforcement, 8x10 cm adhesive fenestration*

Main materials

Absorbent and impervious tri-laminated material (polyethylene + polypropylene)
HSPP absorbent reinforcement (polypropylene + polyethylene)

Pack factor

38 units per carton

Packaging

Paper and film peel pouch
1 polyliner contains all the units in a carton box

Legal Manufacturer

Medline International France SAS
In compliance with EN ISO 13485

Regulation

In compliance with Medical Devices Regulation 2017/745 and UK MDR
MD Class: I sterile
*HP: per EN 13795-1

•Please refer to complete labeling before using the product.

Other Information

NA

Revision: 9
Document Title: FINT.441 - Technical sheet
Document Number: EU1-00004-F-00025

© 2023 Medline Industries LP or one of its subsidiaries. All rights reserved
TDS_PediatricDrape_TB29494CE_EN04

Retention period: End of life and support period



2797 Latex Free

0086