



Medline Industries Ltd
3rd Floor, Quayside Wilderspool
Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

Invisishield™

Probe covers kit

Reference(s)

References	Length (cm)	Width (cm)
ICE5900	122	10
ICE5910	61	15
ICE5920	122	15
ICE5930	244	15

Description

A sterile wrap, 30x30 cm, contains:
- 1 probe cover ⁽¹⁾
- 1 ultrasound gel*, 20g* ⁽²⁾
- 2 rubberbands

Main materials

Clear polyurethane

Pack factor

20 units per carton

Packaging

Paper and film peel pouch
1 polyliner contains all the units in a carton box

Legal Manufacturer

(1) Medline International France SAS.



(2) Parker Laboratories, Inc



Regulation

In Compliance with the Article 22 of Medical Device Regulation (EU) 2017/745
Class of Medical device : PP (Procedure Pack)
* The CE marking of this component is covered by the supplier.

*Please refer to complete labeling before using the product.

Other Requirements

NA

Revision: 9

Document Title: FINT.441 - Technical sheet

Document Number: EU1-00004-F-00025

© 2023 Medline Industries LP or one of its subsidiaries. All rights reserved
TDS_ProbeCovers01_EN04

Retention period: End of life and support period

