

Product Data Sheet



CT Gantry Cover

Ref.:	Product name:	Dimensions:	Weight:	Packaging dimensions:
1231-A	CT Gantry Cover C33	1190 x 1140 x 10 mm	700 g	750 x 450 x 60mm
1271-A	CT Gantry Cover C307	1190 x 1100 x 10 mm	700 g	750 x 450 x 60mm

Intended use	Reusable hygienic CT gantry protection during invasive and non-invasive procedures to keep the CT gantry protected from liquids. Covers the lower half of the CT gantry. May only be used by medical personnel. Non-sterile and not to be used in the sterile field.	
Directives & classification	CE according to MDR – Regulation 2017/745 Class 1 product	
Applied quality assurance	SN EN ISO 13485:2016	
Instructions for use	Enclosed to the product	
Restrictions of use	Do not use the product for purposes other than those described.	
General safety notes	Non sterile. Do not sterilize. Dispose of any product that is cracked or appears damaged Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).	
Materials	Material TPU PA LDPE Cardboard	Component Foil Velcro packaging packaging
Marking	CE	
Storage	Store in dry place / room temperature	

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.