

TOP MASK ULTRA SAFE

DISPOSABLE SURGICAL 3-LAYERS FACE MASK WITH ELASTIC BAND

CONFIGURATION OVERVIEW

- Designed for medical use
- Soft and comfortable
- Water repellent
- Fiberglass-free
- Latex-free
- Good performance of anti micro-dust, anti-spray, anti-virus
- Designed for medical use



QUALITY ASSURANCE

Tests performed: in accordance with EN 14683:2019+AC:2019, EN ISO 11737-1:2018 membrane filtration, iso 22609:2004 splash resistance pressure, EN ISO 14971:2012, EN ISO 15223-1:2016, EN 1041:2008, EN ISO 10993-1:2018/10993-5:2009/10993-10:2013

Conformity to standards: **Medical Device of Class I** according to EU REG MDR 2017/745, Type-IIR, Bacterial Filtration efficiency (BFE) >99%, differential pressure <60

TECHNICAL SPECIFICATIONS

- 3 layers (2 polypropylene layers) + 1 Melt blown layer (middle)
- Outer Layer: 35 g/m², PP, Color
- Middle Layer: 25 g/m², Purpose-made Melt blown, White
- Inner Layer: 25 g/m², PP, nonwoven material, White
- 2 Holding Elastics: urethane elastic fiber
- 1 Nose Barret: polyethylene with iron wire
- Force at break: 6 N

Dimensions mask	Thickness	Weight
17.5 x 9.5 cm	0.6 mm	3.05 g

Item number and colour variation: K01 = 20X50 pcs/cart.

K02 = 40X50 pcs/cart

Blue:	FAMAA25-AKZ-BL-K02	FAMAA24-AKZ-BL-K02/K01	FAMAA52-AKZ-BL-K02
Laguna blue:		FAMAA24-AKZ-LU-K02/K01	
Lilac:	FAMAA25-AKZ-LI-K02	FAMAA24-AKZ-LI-K01	FAMAA52-AKZ-LI-K02
Fresh green:	FAMAA25-AKZ-FG-K02	FAMAA24-AKZ-FG-K02/K01	FAMAA52-AKZ-FG-K02
Fuchsia:	FAMAA25-AKZ-FU-K02	FAMAA24-AKZ-FU-K01	FAMAA52-KZ-FU-K02
Mint green:		FAMAA24-AKZ-EG-K01	
Med blue:	FAMAA25-AKZ-MB-K02	FAMAA24-AKZ-MB-K02/K01	
Pink:	FAMAA25-AKZ-PI-K02	FAMAA24-AKZ-PI-K02/K01	FAMAA52-AKZ-PI-K02
White:		FAMAA24-AKZ-WH-K02/K01	FAMAA52-AKZ-WH-K02
Orange:		FAMAA24-AKZ-HO-K02/K01	
Magic Blue:		FAMAA24-AKZ-MC-K01	
Green:		FAMAA24-AKZ-GR-K02	

PACKAGING

Pieces/dispenser	50	Dispenser/carton	40
Carton size	52x38x34 cm	Carton weight	8.5 Kg
Cartons/layer	4	Layers/pallet	6
Cartons/pallet	24		
Pallet size	80x120x219 cm	Pallet weight	224 Kg

Instructions:

1. Take and face the outer-surface of the mask.
2. While facing the outer-surface of the mask away from you, loop the ear loops comfortably around your ears.
3. Adjust the noseband to fit nicely on the bridge of your nose, conforming to the contour of your face.
4. Now holding the noseband in place, pull the bottom end of the mask under your chin to fully cover your mouth,

your nose, as well as your chin.



Akzenta International SA
Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 91 210 46 00
fax. +41 91 682 85 20

www.akzenta.com
e-mail: business@akzenta.com